

Max Sample P21325

Concept Generation

Function	Concepts									
Wash Breast Pump Parts	Spray with hot water	Spray with hot water and cleaner solution	Leave parts submerged and soak in water solution	Spray with hot water and cleaner solution, spin items during cycle	Soak in hot water and cleaner solution, spin items during cycle	Use scrubbing system (oscillating brushes)	Use scrubbing system, in conjunction with cleaner solution	High water pressure cleaning system	High water pressure cleaning system, with scrubbing system	chemical bath, rinse with water
Dry Breast Pump Parts	Leave in system, air dry	Blow air over parts to decrease drying time	Blow heated air over parts to decrease drying time	Spin parts to dry	Blow air and spin parts simultaneously	Blow heated air and spin parts simultaneously	Expose to sunlight to decrease drying time	Lower pressure of system, water boils off, remove water vapor from chamber	Microwave parts, water boils off, remove water vapor from chamber	Rub parts together vigorously, generate heat and friction